

**Product Name :**  
Demonstration Set Of Bernoulli's Theorem

**Product Code :**  
NLAB-EDUCATIONTVET14009



**Description :**

Demonstration Set Of Bernoulli's Theorem

**Technical Specification :**

The set consists of a Venturi tube made of transparent acrylic material. The pressure is measured at 8 places along the tube with an 8-tube manometric panel.

The flow rate is regulated with a valve. Manometric tube length 300 mm, the tubes can be pressurized with a hand pump. Venturi tube diameters: 25 mm at the inlet and 10 mm at the throat

14 ° upstream cone, 21 ° downstream cone

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003