

Product Name :
Demountable Mechanical Strain Gauge

Product Code :
NLAB -CONCRLAB120029



Description :

Demountable Mechanical Strain Gauge

Technical Specification :

Mechanical strain gauge for the measurement of length variation. Measuring base 100 mm x 5 mm range. Digital gauge, 0.001 mm resolution, output for PC connection (special cable required, see accessories)

Mechanical strain gauge for the measurement of length variation. Measuring base 300 mm x 5 mm range. Digital gauge, 0.001 mm resolution, output for PC connection output for PC (special cable required, see accessories)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003