

Product Name :
Density Basket

Product Code :
NLAB-AGRETLAB160017



Description :

Density Basket

Technical Specification :

Made of Brass / GI with Stainless steel Wire Mesh 6.3mm / 4.75mm size
Ruggedly Constructed, Approximately 20cm dia x 20cm high. Complete with Handle.

Note: Customize sizes also available

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003