

Product Name : Dental Compressor	Product Code : NLABMEDICAL-M3210003
Description : Dental Compressor India, Dental Compressor manufacturer India, suppliers, Medical Dental Compressor India, Hospital Dental Compressor India	
Technical Specification : 1/2 HP Dental compressor seal unit fitted with safety valve for air release excess sorted. Auto cut off system pressure gauge and drain valve air rotor air compressor mounted on high pressure tank entirely powder coated paint Power ½ HP Amperes 4.5 Flow 40 to 50 Liter per minute Speed 1000 rpm Noise 50 to 60 dB (A) Power supply 220V, 50 Hz Tank 25L Pressure 350 MMHG Weight 35 to 37kgs	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003