

Product Name : Dental Double Walled Acrylizer	Product Code : NLABMEDICAL-M32100023
Description : Dental Double Walled Acrylizer India, Dental Double Walled Acrylizer manufacturer India, suppliers, Medical Dental Double Walled Acrylizer India, Hospital Dental Double Walled Acrylizer India	
Technical Specification : Acrylizer Double walled, inside made of stainless steel and out side mild steel sheet painted with attractive stove enamel Temperature is controlled by hydraulic thermostat from 5 c above ambient to 100 C with a sensitivity of + 5 C Acrylizer. Consist Mechanical Timer can be set 1/2 Hour to 2 Hours, brass drain plug, control knob, On/Off switch Acrylizer Inner Size: 300 X 250 X 350 mm	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003