

Product Name :
Dental Double Walled Acrylizer

Product Code :
NLABMEDICAL-M32100023

Description :

Dental Double Walled Acrylizer India, Dental Double Walled Acrylizer manufacturer India, suppliers, Medical Dental Double Walled Acrylizer India, Hospital Dental Double Walled Acrylizer India

Technical Specification :

Acrylizer Double walled, inside made of stainless steel and out side mild steel sheet painted with attractive stove enamel

Temperature is controlled by hydraulic thermostat from 5 c above ambient to 100 C with a sensitivity of + 5 C Acrylizer.

Consist Mechanical Timer can be set 1/2 Hour to 2 Hours, brass drain plug, control knob, On/Off switch

Acrylizer Inner Size: 300 X 250 X 350 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003