

Product Name :
Dental Lathe

Product Code :
NLABMEDICAL-M32100022

Description :

Dental Lathe India, Dental Lathe manufacturer India, suppliers, Medical Dental Lathe India, Hospital Dental Lathe India

Technical Specification :

Power source: 230V~115V or 115V~(on request), 50/60Hz

Weight: 2.4 Kgs.

Dimension: W200 x D140 x H80(mm)

Curing Light Function:

Complete curing restoration with a 75 watt, 12V, halogen lamp

Light guide with 360 degree rotation and can be autoclavable

3 seconds interval to turn on or off the lamp controlled by computer to protect and lengthen use life of lamp

Individual curing time selection with beep tone signal

Automatically turning off fan after cool down period

Provides fuse and thermal cut-off switch to prevent overheating

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003