

Product Name :
Desktop Spectrometer

Product Code :
NLAB-SCHOOLABS21005



Description :

Desktop Spectrometer

Technical Specification :

This tool serves to separate the light into wavelengths makes it possible to identify the components of the light source has

scale diameter. 127mm swivel mounted on base

with fine adjustment screw on telescope and coarse adjusting base by loosening the holder

telescope and collimator equipped with an achromatic objective lens of 32 mm in length. 178mm focus

Ramsden ophthalmic telescope, 15x magnification with crossline pointer, Rack & Pinion focusing system – 6mm slitting slit

mounted on the telescope. and collimeter

Includes prism holder, grating holder, goggles and tommy -bar for adjusting the axle

all housed in a wooden case

features the same as PL0291, but with a 150mm dial diameter scale mounted on the base.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003