

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |                                               |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------|
| <b>Product Name :</b><br>Diagnostic Audiometer                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | <b>Product Code :</b><br>NLABMEDICAL-M3220013 |
| <b>Description :</b><br>Diagnostic Audiometer India, Diagnostic Audiometer manufacturer India, suppliers, Medical Diagnostic Audiometer India, Hospital Diagnostic Audiometer India                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |                                               |
| <b>Technical Specification :</b><br>Input: Tone, speech-built in & Ext. Mike<br><br>Output: Left, Right, Bone<br><br>Freq: Air (-20 db to 120 db) 250 Hz to 8 Khz - 5 db Step. Bone (-20 db to 80db) 250 Hz to 6 Khz- 5 db Step.<br><br>Masking: Wide Band-10db Step., 100 db Max.<br><br>Tests: Pure Tone, Speech, Tone Decay (counter)<br><br>Tone: Continuous & Pulse<br><br>Screen: L.C.D.(with contrast control)<br><br>Head Phone: DR 59, TDH 39/49/Optional<br><br>Bone: B-71 Radio ear<br><br>Indication: V.U. Meter for Speech<br><br>Voltage: AC 230 / 50 Hz, DC 12V (IA) Operation<br><br>Weight: 3.5Kg. (With Shoulder Bag)<br><br>Dimension: 25 x 21 x 7.5 cm |                                               |

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com  
**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

## Naugralabequipments