

<b>Product Name :</b> Diagnostic Insulation Resistance Tester 15 KV	<b>Product Code :</b> NLAB-TESTINGLABM37002
--	--

--

**Description :**

Diagnostic Insulation Resistance Tester 15 KV

**Technical Specification :**

Description

Diagnostic tests for polarization index, dielectric absorption ratio, step voltage, dielectric discharge and ramp testing.

For the diagnostic testing and maintenance of high voltage electrical equipment

Measures insulation resistance to 30 T?.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003