

Product Name :
Dial Gauge

Product Code :
NLAB-LBEMALAB230028



Description :

Dial Gauge

Technical Specification :

Attractive Shape And Elegant Design Complimented By Uncompromising Performance !

With 0.01mm Reading ,Range 10mm & Graduation 0-100

Well-Balanced Frames And Handgrips Coupled With Lightweight Make These Gauges Easy To Use For The Measurement.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003