

Product Name :
Dial Thermometer

Product Code :
NLAB-LBEMALAB230023



Description :

Dial Thermometer

Technical Specification :

DIAL THERMOMETERS work on principle of bi metal expansion. Scale or range of temperature which DIAL THERMOMETER can measure may vary from -200 0 C to 650 0 Centigrade. DIAL THERMOMETER consists of a sharp metal indicator which is fixed on a circular temperature scale. There are two main components , one is a needle and another is a temperature sensor. Few more types of DIAL THERMOMETERS are Liquid type, Gas filled and Va pour tension dial thermometers.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003