

Product Name :
Differential Voltage Sensor

Product Code :
NLAB-VOCATIONALTVET28008



Description :

Differential Voltage Sensor

Technical Specification :

Supply voltage: 5 VDC

Input Impedance (to ground): 10 M ?

Input voltage range: ± 6.0 V

Maximum voltage on any input: ± 10 V

Differential Impedance: >20 M ?

Resolution: 3.1 mV

Linearity: 0.01%

We are well-known manufacturers, OEM suppliers of Differential Voltage Sensor for Data Acquisition System. Contact us for high quality Differential Voltage Sensor for Data Acquisition System for schools lab, college lab, universities, research labs, various teaching and workshop training laboratories and industries in India.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003