

Product Name :
Physics Lab Digimatic Caliper

Product Code :
NLAB-EDUSCHOOLEQ158007



Description :

Physics Lab Digimatic Caliper

Technical Specification :

Hardened and tempered stainless steel digital calliper with satin finish

Large LCD display with 11mm high figures for ease of reading

Material: Stainless Steel

Measuring range: 0-300mm

Can be measured OD, ID, depth and steps, Would have a thumb roller

Buttons: On/Off, Zero

Display: LCD

Power Supply: Battery

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003