


**Product Name :**  
Digital Bulk Density Apparatus

**Product Code :**  
NLAB-SCIENTIFICEQ28007



**Description :**

Digital Bulk Density Apparatus

**Technical Specification :**

4-digit microprocessor based Pre-set counter for the setting of exact no of strokes.

A common holder for both the cylinder drops down.

Easy to lock cushioned holders are provided to hold the cylinder.

In power failure, the non-volatile memory of the microprocessor remembers the last mode of operation and complete the test after the power supply is resumed.

Technical Specification:

Cylinder Capacity: 50 Ml. x 2 Nos.

of Drops : Adjustable between 1 to 9999 drops.

Display : Digital 12 mm Red LED

Drop Height : 25 mm approx.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003