

Product Name :
Digital Bulk Density Apparatus

Product Code :
NLAB-SCIENTIFICEQ28007



Description :

Digital Bulk Density Apparatus

Technical Specification :

4-digit microprocessor based Pre-set counter for the setting of exact no of strokes.

A common holder for both the cylinder drops down.

Easy to lock cushioned holders are provided to hold the cylinder.

In power failure, the non-volatile memory of the microprocessor remembers the last mode of operation and complete the test after the power supply is resumed.

Technical Specification:

Cylinder Capacity: 50 ml. x 2 Nos.

of Drops : Adjustable between 1 to 9999 drops.

Display : Digital 12 mm Red LED

Drop Height : 25 mm approx.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003