

Product Name :
Digital Clamp Meters Machines

Product Code :
NLAB-TESTINGLABM18001



Description :

Digital Clamp Meters Machines

Technical Specification :

Can measure AC and DC currents up to 2000A.

- Can measure the Voltage and Current in both very low and high power circuits.
- Accurate True-RMS reading of AC current or voltage with distorted waveform.
- Output terminal for connection to recorders.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003