

**Product Name :**  
Digital Clampmeter 300A AC

**Product Code :**  
NLAB-TESTINGLABM18003



**Description :**

Digital Clampmeter 300A AC

**Technical Specification :**

Maximum Conductor Size: 32 millimeter

Large LCD for easy reading, data hold switch freezes reading

Low battery indication, continuity buzzer

Tough ABS plastic housing, diode check

Accessories: Test leads, carrying case, manual, drop proof wrist strap

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003