

Product Name :
Digital Clampmeter

Product Code :
NLAB-TESTINGLABM18004



Description :

Digital Clampmeter

Technical Specification :

Digital Clampmeter:-

3 3/4 Digital Clampmeter 1000A – AC

Features:

Rish clamp 300A is a highly innovative design for features those increases safety and comfort of user. Rotating clamp jaws facilitate the measurement at physically awkward positions, vertical bus bars, conductors placed at positions difficultla access

Clamp jaws can be opened or closed with the trigger placed at bottom side away from the jaws This allows the user to place his/her hand at safer distance from live conductor. This greatly reduces exposure of human beings to electrical shocks

Location and design of trigger eliminates fatigues caused by single finger operation. It allows spreading the force required to open the jaws over more Than one finger to ensure comfortable operation

Comfortable operation of push buttons and function selector switch, in adverse field conditions

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003