

**Product Name :**  
Digital Colposcope

**Product Code :**  
NLABMEDICAL-M3470001



**Description :**

Digital Colposcope India, Digital Colposcope manufacturer India, suppliers, Medical Digital Colposcope India, Hospital Digital Colposcope India

**Technical Specification :**

High-person optics

Reliance for in-office surgical application

Also excellent for operation theater and lab work

Very economically priced

Also conveniently moveable for one place to another, even by lightly-built/lady gynaecologists

Gynecological (Including and uterine) diagnosis

Detection of uterine cancer, or uterine erosion, tumors, congestion of inflammation.

Very accurate Biopsy of uterus or

Specifications:

Tilting Adjustment: 30 Upward, 30 Downward

Working Range(Fine Focusing: 50mm

Lamp: 15V/150W,24V/250W

Total all-inclusive weight: 22kg.

filter: Red Free Green

Individual Components:

Eye Piece: 15 x,12.5x

Objective Lens: 200mm,300mm,400mm

Field of vision: 32mm

Magnification: 5x,10x,20x-3 step Model maxcp 1 5x7.5x,12.5x,30x-5 Step Model-VICS-30 5x to 25x continuous Zoom Model-maxcp-1

Cross Slide: 50mm

Light Source: Halogen Fibre, Optic, Single Or Double outlet

Optional Accessories: Beam Splitter With CCD Camera Attachment

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003