

Product Name :
Digital Conductivity Meter

Product Code :
NLAB-ENVRTLAB200005



Description :

Digital Conductivity Meter

Technical Specification :

Display 0.5 Digit LED.

Range 0 to 1000 Mhos / cm.

Accuracy ± 1 digit.

Temp. Compensation 0 to 50 deg.C.

Adjustable

Sensor Conductivity Cell.

Power 10% 50 Hz. AC

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003