

Product Name :
Digital Frequency Counter Meter

Product Code :
NLAB-TESTINGLABM21002



Description :

Digital Frequency Counter Meter

Technical Specification :

Technical Specifications

Measuring range : 10MHz, 100MHz, 2500MHz (3 ranges)

Accuracy : $\pm(4 \text{ PPM} + 1\text{d})$

Time base circuit : 4.194Mhz quartz crystal

Display : 0.5" LCD and 8 digits

Intelligent Functions : Frequency, Period, Multi-resolution, Data hold, Relative measurement Data record (Max. Min. Average reading)

Input connector : BNC, female

Case : Durable and strong ABS plastic housing

High sensitivity for VHF and UHF frequency measurement

Auto power off and manual power off function

Power : DC 6 volt battery

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003