

<b>Product Name :</b> Digital IC Tester Trainer	<b>Product Code :</b> NLAB-EDUSCHOOLEQ197008
<b>Description :</b>	
Digital IC Tester Trainer	
<b>Technical Specification :</b>	
Study and Construction of Various Digital Electronics Lab Experiment	
TECHNICAL SPECIFICATIONS:	
Power Supplies:	
Dc Power Supply IC Regulated +5V DC, 1A.	
Operated On Mains Power 230v, 50hz +10%	
20 Pin Zip Socket (3Nos.)	
40 Pin Zip Socket (1Nos.)	
Logic Inputs & Outputs 12 Bit	
4 Bit BCD To 7 Segment Indicators (3Nos.)	
Logic Low Frequency Mono Pulsar.	
Low Frequency Clock Generator.	
Buzzer as Continuity Tester.	
Variable Resistors (2Nos.)	
Resistors and Capacitors.	
SALIENT FEATURES:	
Front Panel Built With High Class Insulated Printed Circuit Board Sheet With Well Printed Circuits And Symbols.	

Fuse For Short Circuit Protection

Instruction Manual.

Connections are brought out through 2mm Brass Sockets.

Patch Cords 2mm.

The trainer is housed in ABS Plastic cabinet.

Size of the trainer set 12"x10"

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003