

**Product Name :**  
Digital IC Trainer

**Product Code :**  
NLAB-ENGINEERINGLB13002



**Description :**

Digital IC Trainer

**Technical Specification :**

Digital IC Trainer Lab Equipments Manufacturer, Suppliers & Exporters

Digital IC Trainer (TTL) has been designed for verification of the truth table of logic gates and address etc.

Students can be trained about the digital ICs by using this board. The fundamentals of the digital ICs can easily be understood by performing required experiments for verification of truth table of logic gates. Quite large area of the bread board has been provided on the panel for ICs . Students can make the circuits easily on the bread board with the help of other accessories which have been provided on the panel. Practical experience on this board carries great educative value for Science and Engineering Students

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003