

**Product Name :**  
Digital Impact Test Machine

**Product Code :**  
NLAB-TESTINGLABM40003



**Description :**

Digital Impact Test Machine

**Technical Specification :**

The Digital Impact Testing Machine mainly used to determine the anti-impact capability and strengthen impact toughness of non metallic materials such as nylon, glass fiber reinforced plastic, ceramic & cast stone, electrical insulating materials etc.

Sample and support are easy centering

The pendulum of the machine can be raised or released automatically by button pressed

Digital LCD screen control and value display

The machine comply with GB / T1843, ISO 180, GB / T21189, GB / T2611 and other international standard

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

