

**Product Name :**  
Digital Iso Block Dry Bath

**Product Code :**  
NLAB-EDUSCHOOLEQ45002



**Description :**

Digital Iso Block Dry Bath

**Technical Specification :**

Iso Block dry bath is ideal for incubation of samples at two different temperatures, in a single compact unit.

The hinged (and removable) lid is also designed with a dividing wall.

The two chambers are separated by an insulated wall and are controlled independently, allowing IsoBlock to serve as two dry baths in one.

When closed, this dividing wall acts as a barrier, preventing air transfer between chambers.

In addition, the unique Quick-Flip universal blocks accept tubes from 0.5 to 2.0ml and, with just a quick flip to the other side, they also accommodate 0.2ml tubes or PCR strips.

The timer function is non-mandatory; Iso Block can be used for continuous heating.

The control panel features large digital displays of both time and temperature.

Both settings of a block can be adjusted at any point, without affecting the settings of the other block.

Various blocks are available for accommodating nearly all popular tube sizes

Specification:-

Ambient Temperature Range	
---------------------------	--

Temperature Accuracy:	
-----------------------	--

Temperature Uniformity:	
-------------------------	--

Temp. Controller Increments:	
------------------------------	--

Offer: 1 min. to 99 hr 59 min.	
--------------------------------	--

Dimensions 6.3 x 14.5 x 5 in. / 16 x 36.5 x 13 cm	
--	--

Weight: 4.5kg	
---------------	--

Power: 230V, 50-60Hz	
----------------------	--

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003