

Product Name :
Digital Lab Station

Product Code :
NLAB-VOCATIONALTVET14004



Description :

Digital Lab Station

Technical Specification :

Digital Lab Station The board consists of the following built -in parts. +5V D.C. at 200mA,
IC Regulated Power Supply internally connected. Nine,
3-input AND gates whose outputs is connected to an inverter (NOT gate) to give 3-input NAND gates.
A clock generator with a repetition frequency of 5 Hz.
Two lamps (LEDs) driver circuits each of which individually drivers a lamp (LED).
Switches ofr logic selection. Lamps (LEDs) for visual indication of status.

Electrical Work Bench Lab Digital Lab Station Manufacturer,Supplier and Digital Lab Station Exporter in India

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003