

**Product Name :**  
Digital Lab Station

**Product Code :**  
NLAB-VOCATIONALTVET14004



**Description :**

Digital Lab Station

**Technical Specification :**

Digital Lab Station The board consists of the following built -in parts. +5V D.C. at 200mA,  
IC Regulated Power Supply internally connected. Nine,  
3-input AND gates whose outputs is connected to an inverter (NOT gate) to give 3-input NAND gates.  
A clock generator with a repetition frequency of 5 Hz.  
Two lamps (LEDs) driver circuits each of which individually drivers a lamp (LED).  
Switches ofr logic selection. Lamps (LEDs) for visual indication of status.

Electrical Work Bench Lab Digital Lab Station Manufacturer,Supplier and Digital Lab Station Exporter in India

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003