

Product Name :
Digital Load Frame

Product Code :
NLAB-TESTINGLABM109055



Description :

Digital Load Frame

Technical Specification :

It is an advanced version of mechanical load frames. In this load frame, rate of strain is controlled through combination of electronics and micro stepper motor. The unit operates on 220Volts 50Hz single phase supply.

Specification:

- Capacity of the load frame 50kg
- Vertical Clearance 180mm
- Horizontal Clearance 135mm
- Platen Diameter 150mm
- Platen Speed 99mm/min
- Weight App. 160kg

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003