

Product Name : Digital Microscope Pro	Product Code : NLAB-EDUSCHOOLEQ173080
Description :	
Digital Microscope Pro	
Technical Specification :	
<p>The Digital Microscope Pro is an easy to use, low-power microscope with a 5.0 MP sensor for capturing photos and video of your discoveries.</p> <p>At low magnification, the Handheld Digital Microscope Pro is ideal for viewing stamps, coins, bugs, plants, rocks, skin, gems, circuit boards, and more. You can even view traditional slides at high magnification.</p> <p>Specifications:-</p> <p>Professional, adjustable height stand included for hands-free operation.</p> <p>4 foot USB 2.0 cable for easy maneuverability when viewing large surfaces.</p> <p>5-Element IR cut high-quality glass lens ensures sharper images.</p> <p>20x to 200x powers, great for low-power observation of 3D specimens.</p> <p>Windows and Mac compatible.</p> <p>Intuitive software with measuring features.</p> <p>True 5MP sensor to capture and save high-resolution images and videos of your specimens.</p>	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

