

**Product Name :**  
Digital Multimeter

**Product Code :**  
NLAB-ELECTRONICSAB190001



**Description :**

Digital Multimeter

**Technical Specification :**

Digital Multimeter We are Leading Supplier of Digital Multimeter in India. We also bulk supply Digital Multimeter for laboratory tenses. Features 3 1/2 Digit Manual Ranging Temperature Measurement Capacitance Measurement Data hold and auto power off function Technical Specifications : DC Voltage : Range : 200mV ~ 600V Accuracy :  $\pm(0.5\% + 1\text{digit})$  AC Voltage : Range : 2V ~ 600V Accuracy :  $\pm(0.8\% + 3\text{digits})$  DC Current : Range : 2mA ~ 10A Accuracy :  $\pm(0.8\% + 1\text{digit})$  AC Current : Range : 20mA ~ 10A Accuracy :  $\pm(1.2\% + 5\text{digits})$  Resistance : Range : 200 ~ 20M $\Omega$  Accuracy :  $\pm(1.0\% + 2\text{digit})$  Diode test open circuit : 2.8 volt Capacitance : Range : 2nF ~ 200 $\mu$ F Accuracy :  $\pm(4.0\% + 3\text{digits})$  Temperature : Range : -40°C ~ 1000°C Accuracy :  $\pm(1\% + 7\text{digits})$  Power : DC 9 V battery

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003