

Product Name :
Digital Multimeter Machines

Product Code :
NLAB-TESTINGLABM19001



Description :

Digital Multimeter Machines

Technical Specification :

This DT92 Advanced Series Digital Multimeter is a compact precision, battery operated LCD display 3-1/2 or 4-1/2 digits Digital Instrument.

High Accuracy,

Large rotational folded LCD display,

button release lock,

Digital Height 25mm,

Single 32 position rotary switch for FUNCTION and RANGE selection,

Allows fast and convenient operation.

Cuvilinear mode soft case.

Coloured indication jack with fully protection test leads.

Lower average power Auto-Power Off

Data hold for easy reading

DC Voltage

AC Voltage

DC Current

AC Current

Transistor test

Capacitance Test

OHM

Diode Test

Continuity Buzzer

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003