

Product Name :
Digital Multimeter

Product Code :
EC154



Description :

Digital Multimeter

Technical Specification :

Digital Multimeter Auto-Ranging Digital Multimeter 3½ Digit, 2000 count LCD display with backlight AC/ DC voltage, AC/DC current, resistance and frequency measurement Diode and continuity test functions Weight: 8 ounces Dimensions: 2 x 10 x 5 inches Batteries: 1 AAA batteries required. (included)

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003