

Product Name : Digital Multipara Monitors	Product Code : NLABMEDICAL-M3360003
Description : Digital Multipara Monitors	
Technical Specification : 8.4 = color TFT display. Suitable for adult, pediatric and neonatal patient. Basic parameters (ECG, SpO2, NIBP) in a durable case for bedside monitoring and transport monitoring. Available and visual alarms with adjustable alarm ranges. Net-workable data management and storage capacity.	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com
Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003