

Product Name :
Digital Nephelometer

Product Code :
NLAB-LABORATORYEM31001



Description :

Digital Nephelometer

Technical Specification :

Range: 0 to 200 JTU 0 to 1000 JTU 2 Ranges

Resolution: 1 JTU

Accuracy: +3% FS + 1 digit

Sample System: 30 mm clear glass tube

Light Source: 6.8 V, Tungsten Lamp

Detector: Photodiode

Display: 3 1/2 digit LED display

Power: 220V+ 10% AC 50 Hz

Size (LxBxH): 305 x 305 x 165 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003