

**Product Name :**  
Digital PH, Conductivity and Temperature Meter

**Product Code :**  
NLAB-LABORATORYEM15002



**Description :**

Digital PH, Conductivity and Temperature Meter

**Technical Specification :**

Digital pH, Conductivity & Temperature Meter

Highly Stable and Accurate

3½Digit Red LED Display

Cell Constant Adjustment Facility

Rugged & Reliable Design

Digital pH, Conductivity & Temperature Meter is an elegant rugged and most economical instrument for fast and accurate pH, Conductivity and Temperature measurements in any laboratory.

This instrument is useful for monitoring acidity & alkalinity, conductivity & temperature of natural water, sea water, drinking water, treated water, waste water, brine solution, soil and other chemical solutions.

3½ digit, seven segment bright red LED display & solid state IC circuitry makes it versatile and reliable.

It operates on 230V AC, 50 HZ.

These are extremely useful instruments for chemical laboratories, pollution and environmental testing, public health engineering, sugar cement, paper industries, pharmaceuticals & boiler water analysis, agriculture and soil labs, brewery and transport under takings, fertilizer plants and petroleum refineries etc.

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003