

Product Name :
Digital PH Meter

Product Code :
NLAB-ENVRTLAB200001



Description :

Digital PH Meter

Technical Specification :

Display 3 Digit LED.

Range ± 14.00, mV 0 to ± 1999.

Resolution 0.1 mV.

Accuracy ± 0.1, mV ± 1 mV.

Temperature Compensation, Manual 0 to 100 deg. C.

Storage 120% rol

Power 10% 50 Hz AC.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003