

**Product Name :**  
Digital Point Load Test Apparatus

**Product Code :**  
NLAB-ROCKTLAB170009



**Description :**

Digital Point Load Test Apparatus

**Technical Specification :**

Digital Point Load Test Apparatus consists of a 60 kN capacity load frame with a hydraulic loading ram driven by a hand pump. The frame is adjustable for testing of samples up to 102mm diameter. A ruler assembled on the frame allows the direct measurement of the distance between the conical platens before and after the test. The compression load is measured by a pressure transducer connected to an advanced digital display unit assuring the best accuracy and resistance to the failure shocks.

Digital Point Load Index Tester

**Features:**

Load range : 0 – 60 kn

Digital display : 2 x 16 characters

Resolution : 32.000 div.

Accuracy :  $\pm 1\%$

Load pacer included

Load measurement in both kn and mpa

Serial port for pc connection

## Naugralabequipments

**Website:** [www.naugralabequipments.com](http://www.naugralabequipments.com), **Email:** [sales@naugralabequipments.com](mailto:sales@naugralabequipments.com)

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003