

Product Name :
Digital Point Load Tester

Product Code :
NLAB-VOCATIONALTVET17012



Description :

Digital Point Load Tester

Technical Specification :

The applied load is measured by a high precision electric load cell with a digital display unit range 0-56kN proving:-

65.000 divisions

0.001 kN resolution

Linearity: 0.05%

Hysteresis: 0.03%

Repeatability: 0.02%

The strength index (IS) is got by the formula $P : D^2$ where P is the strength and D the space between the two conical points.

Supplied complete with wooden carrying case, goggles, accessories.

Dimensions: 400 x 530 x 720 mm

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003