

Product Name :
Digital Potentiometer

Product Code :
NLAB-ELECTRONICSAB190005



Description :

Digital Potentiometer

Technical Specification :

Digital Potentiometer is a precision instrument for potentiometric measurements. Results are displayed in milli volts on a direct digital readout by 3½ digit LED display. It is supplied with a set of 4 electrodes viz. Reference, Platinum, Silver & Glass. It is extremely useful potentiometric, ORP and acid base titrations.

FEATURES :

Highly Stable & Accurate
4 Electrodes for different Titrations
3½ Digit Led Display
Rugged & Economical.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003