

<b>Product Name :</b> Naugra Lab Digital Respinometer	<b>Product Code :</b> NLAB-LABORATORYEM21013
--	---

--

**Description :**

Naugra Lab Digital Respinometer

**Technical Specification :**

Polished Stainless Steel tank, with capacity of 7 manometers a side.

Digital temperature display with resolution of 0.1 deg. C.

Temperature control to within +/-0.5 deg. C.

Stainless Steel immersion heater.

Stirrer for uniform temperature.

Manometers mounted on graduated scale for reading pressure changes.

Flask with center well & side flask attached to manometer.

Oscillation of manometer approximately 90 rpm.

Powder coated exterior, with wheels for easy mobility.

Capacity 7 manometers, temperature range of 5 deg.

C above ambient to 60 deg. C

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003