

**Product Name :**  
Digital Sound Wave Kit with Resonance Tubes

**Product Code :**  
NLAB-SCHOOLABS15005



**Description :**

Digital Sound Wave Kit with Resonance Tubes

**Technical Specification :**

To study the resonance or resonance of sound in the pipe ends. and the relationship of frequency to the number of standing waves produced when the medium is air. As well as in trials for the speed of sound in air

containing

Signal Generator. Generates a sinusoidal sound wave in the frequency range 1 kHz – 3kHz. The volume can be adjusted. digital display screen There is a connection for 2 speakers.

Resonance tube, total length of the device approx. 105 cm., showing a distance scale of 80 cm., with a sliding bow. One end is attached to the speaker. Resonance tube Ø 60 mm. The speaker is equipped with Bluetooth (Bluetooth) can be connected to the Function Generator application on a smartphone or other device.

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003

