

Product Name :
Digital TDS Meter

Product Code :
NLAB-ENVRTLAB200004



Description :

Digital TDS Meter

Technical Specification :

Display LED.

Power 10% 50 Hz.

Range 1000 ppt in 5 Ranges.

Resolution

Accuracy ± 1 Digit.

Sample Size

Eligibility Type.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003