

Product Name :
Digital Thermometer

Product Code :
NLAB-EDUCATIONTVET15029



Description :

Digital Thermometer

Technical Specification :

For temperature measurement and temperature difference with two NiCr-Ni thermocouples.

Liquid crystal display for displaying two temperatures. T1, T2 or T1-T2 display

Selectable ° C or ° F units

Storage of MIN & MAX extreme values

HOLD function

Temperature:

Measuring range: -50 to +1,350 ° C

Resolution: 0.1 ° C up to 200 ° C, 1 ° C outside

Accuracy: +/- 0.2% full scale

Probes: NiCr-Ni thermocouple (max. 2 thermocouples)

Special thermocouple connector

NiCr-Ni thermocouples not included:

Submersible up to 1000 ° C; STAINLESS STEEL

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003