

**Product Name :**  
Digital to Analog Converter Module

**Product Code :**  
NLAB-VOCATIONALTVET14005



**Description :**

Digital to Analog Converter Module

**Technical Specification :**

Software-programmable ranges 12-bit resolution Connect channels in series for output voltages up  $\pm 60$  V per module 250 Vrms working isolation per channel NI-DAQmx measurement services software to simplify configuration and measurement

Electrical Work Bench Lab Digital to Analog Converter Module Manufacturer, Supplier and Digital to Analog Converter Module Exporter in India

## Naugralabequipments

**Website:** www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

**Address:** 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003