

Product Name : Digital Triaxial Accessores	Product Code : NLAB-TESTINGLABM109061
Description :	
Digital Triaxial Accessores	
Technical Specification :	
1 Lateral extensometer (submersible type).It is a ring type extensometer which exerts a bare nominal pressure on the specimen. Extension/contraction is measured by submersible L.V.D.T and the parameter is directly displayed on the monitor in engineering units. *	

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg,Ambala Cantt,Haryana,India. **Phone:** +91-9896600003