

Product Name :
Digital Turbidity Meter

Product Code :
NLAB-ENVRTLAB200003



Description :

Digital Turbidity Meter

Technical Specification :

Display LED.

Power ± 10% 50 Hz.

Range 1000 NTU in 2 ranges.

Resolution

Accuracy ± 1 Digit.

Sample System Glass Test Tubes.

Light Source Tungsten Lamp.

Detector Photodiode

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003