

Product Name :
Digital Vernier Caliper

Product Code :
NLAB-LBEMALAB230031



Description :

Digital Vernier Caliper

Technical Specification :

Hardened stainless steel construction with 4 way measurement: outer diameter, inner diameter, depth, and step.

Readout switches between inches and metric measurements with a push of the Inch/Metric button.

The Zero/ABS button allows for the display of the absolute coordinate or true position from the origin point.

Supplied in fitted plastic case.

package include 1 x Mitutoyo 500-196 Micrometer

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003