

Product Name :
Direct Reading Spectrometer

Product Code :
NLAB-SCHOOLABS21008



Description :

Direct Reading Spectrometer

Technical Specification :

Can spectra of compounds in the wavelength range 400 - 700 nm reading the spectrum of light directly and see. Sharp spectral lines with good quality rating of 600 lines/mm. Slots for cuvettes. Slit length can be adjusted. The instrument provides spectral samples of 16 common elemental wavelengths,

including a spectrometer, 5 transparent test tubes.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003