

Product Name :
Direct Shear (Motorized 1 Speed)

Product Code :
NLAB-TESTINGLABM109047



Description :

Direct Shear (Motorized 1 Speed)

Technical Specification :

Same as above but electrically operated to give a single constant rate of stain of 1.25mm/min. Suitable for 230V A.C. Weight 175kg. approx.

Naugralabequipments

Website: www.naugralabequipments.com, **Email:** sales@naugralabequipments.com

Address: 6148/6, Guru Nanak Marg, Ambala Cantt, Haryana, India. **Phone:** +91-9896600003